



Chloride Titrant, Silver Nitrate 0.0493 N (100 mL)

Продукт №:

2349832

Опасен

BGN Цена:

Свържете се с нас

Налично



Refer to the Test Kit Replacement Reagent Section for proper kit replacements Silver nitrate, 0.0493 N. Measured dropping bottle.

For determination of chloride by drop-count titration. Silver nitrate, 0.0493 N. 100 mL Marked Dropping Bottle.

Спецификации

Метод, име:

Silver Nitrate Titration

Обем / опаковка, размер:

100 ml

Платформа:

Разтвор