



Calibration kit, Stablcal turbidity standards, 2100A turbidimeter, sealed vials

Продукт №:

2659105

Опасен

BGN Цена:

Свържете се с нас

Налично



Stablcal Standards Calibration Kit, for 2100A Laboratory Turbidimeter. USEPA accepted for calibration purposes. Stablcal Stabilised Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. Kit includes Sealed Vials - 1.0, 10, 100 and 1000 NTU

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Спецификации

Концентрация: 1,0, 10, 100, 1000 NTU

Обем / опаковка, размер: Set of 4