



Calibration kit, Stablcal turbidity standards, 2100AN / AN IS turbidimeter, 100 mL bottles

Продукт №:

2659510

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Ready for use, stabilised Formazine Primary Standards in bottles, Set: <0.1, 20, 200, 1000, 4000, 7500 NTU

Спецификации

Инструмент:	2100 AN / AN IS
Концентрация:	< 0,1, 20, 200, 1000, 4000, 7500 NTU
Обем / опаковка, размер:	Set of 6
Описание:	Stablcal Calibration Kit for 2100 AN / AN IS (100 mL Bottles)
Параметър:	Turbidity
Платформа:	Разтвор
Условия на съхранение:	5 °C - 25 °C