



Stablcal Ampule Calibration Kit, Formazin, for 2100Q/QIS portable turbidimeter

Продукт №:

2971205

Опасен

BGN Цена:

Свържете се с нас

Налично



Calibration Standard kit for 2100Q/QIS portable turbidimeter. Ampule kit contains Formazin standards in concentrations: 10, 20, 100, 800 NTU.

Precise, Pre-mixed Formazin Solutions

Based on Formazin Primary Standards

Enhanced Stability Saves Time and Labor

Confidence in results

Спецификации

Концентрация: 20, 100, 800, 10 NTU (for verification)