



Secondary Turbidity Standards Kit for 2100 Portable Turbidimeters

Продукт №:

2464105

Опасен

BGN Цена:

Свържете се с нас

Налично



Kit includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU.

Gelex Secondary Standards Kit, for 2100P and 2100Q Portable Turbidimeters (includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU).

Спецификации

Обем / опаковка, размер:

Set of 3