



Fluoride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Продукт №:

258999

Опасен

BGN Цена:

Свържете се с нас

Налично



Ionic Strength Adjustor (ISA) for Fluoride determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Fluoride level readings. Eliminate most interferences. Pack of 100 powder pillows.

Спецификации

EPA compliant:	Да
Package settings:	Powder Pillows
Измерване, диапазон:	0.01 mg/L to saturated
Концентрация:	ISA
Обем / опаковка, размер:	100 /pk
Параметър:	Fluoride