



Amino Acid F Reagent for Silica analysis, 2.9 L

Продукт №:

2353103

Опасен

BGN Цена:

Свържете се с нас

Налично



Amino Acid F Reagent (2.7 L) used in the Series 5000 Silica Analyser.

Amino Acid F reagent is supplied in two parts to ensure a long shelf life: a 3-liter bottle of Amino Acid F Solution and a bottle of Amino Acid F Powder. Powder must be added to the solution immediately before installation in the analyser.

Спецификации

Измерване, диапазон:

N/A

Инструмент:

5500sc Silica

Обем / опаковка, размер:

1 x 2.7 L & 1 x 308 g

Параметър:

Silica

Условия на съхранение:

Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.