



APA6000 Alkalinity Cleaning Solution

Продукт №:	2697053
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

The APA6000 Alkalinity Process analyzer uses five reagents each of which is supplied in a 1L bottle designed to provide reagents for 30 days of continuous use.

APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample.

APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints of the titration.

APA Alkalinity Standards 1 and 2 (6001100) are used together to perform a multiple point calibration.

The fifth reagent, APA Alkalinity Cleaning Solution (269703), cleans the system during the Prime and Instrument Clean functions.

Спецификации

Измерване, диапазон:	N/A
Инструмент:	APA6000
Обем / опаковка, размер:	1 L
Параметър:	Alkalinity
Условия на съхранение:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.