



CDC566T General-purpose conductivity cell, 4-pole

Продукт №: E61M010
BGN Цена: Свържете се с нас
Не е показана дата на доставка

4-pole conductivity electrode with MAB6 plug and epoxy body. Cell constant 1cm^{-1} . With temperature sensor. Max. temperature 80°C, min. immersion depth 35mm.

Спецификации

Conductivity Cell:	CDC566T
Гаранция:	12 месеца
Електрод:	4-pole, non-platinised
Специална функция:	Conductivity electrode with MAB6 plug and epoxy body.