



pHC2011-8 Combined pH electrode, alkal. Samples, Red Rod, BNC

Продукт №:

E16M317

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Long-term accuracy and reproducibility with Red Rod technology

pH electrode for alkaline samples. With Red Rod reference system. Measuring range 0-14 pH, 0-100 °C. Liquid junction: porous pin. With BNC plug.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Спецификации

Гаранция:	12 месеца
Датчик, тип:	Glass with membrane for reduced alkaline (sodium) error
Дължина:	103 mm
Дължина на кабела:	1 m
Електрод, тип:	pHC2011-8
Измерване, диапазон:	0 - 14 pH
Материал:	Glass
Разтвор за пълнене:	Refillable: Saturated KCl
Специална функция:	Special glass membrane for alkaline samples and Red-Rod reference system
Съединение на сондата:	Porous pin
Температурен диапазон:	0 - 100 °C
Термистор:	None (external required)
Тип на сондата:	Laboratory combination pH electrode with Red-Rod reference system

Задължителни аксесоари

- Разтвор за пълнене, референтен, наситен KCl (KCl.L), 100 mL (Item S21M002)
- Разтвор за пълнене, референтен, кристали KCl (KCl.C), 15 g (Item S21M001)