



PHC2085 Combined pH electrode, Red Rod, Temp. sensor

Продукт №:

E16M501

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: porous pin. With temperature sensor.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Спецификации

Гаранция:	12 месеца
Датчик, тип:	Glass
Дължина:	103 mm
Дължина на кабела:	1 m
Измерване, диапазон:	0 - 12 pH
Материал:	Sensor body: Epoxy
Размери (диаметър x дължина):	12 mm x 103 mm
Разтвор за пълнене:	Refillable: Sat. KCl
Специална функция:	Red-Rod reference system
Съединение на сондата:	Porous pin
Температурен диапазон:	-10 - 100 °C
Термистор:	NTC 10 kohm
Тип на сондата:	Laboratory