



Water, deionised and demineralised, 100mL

Продукт №:	27242
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S}/\text{cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g}/\text{L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g}/\text{L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Спецификации

Платформа:	Разтвор
------------	---------