



Water, deionised and demineralised, 500mL

Продукт №:	27249
BGN Цена:	Свържете се с нас
Налично	

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S/cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g/L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g/L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Спецификации

Обем:	500 mL
Описание:	Water, Deionized
Платформа:	Разтвор