



Chlorophosphonazo Indicator Solution Pillows, pk/100

Продукт №:

2589599

Опасен

BGN Цена:

Свържете се с нас

Налично



For the determination of Ultra Low Range Hardness, 8-1000 µg/L CaCO₃. Chlorophosphonazo indicator solution is packaged in unit-dose solution pillows.

For ultra low range hardness determination by the Chlorophosphonazo method. Hach Method 8374. Range: 4 - 1000 µg/L Ca & Mg as CaCO₃. Chlorophosphonazo solution reagent packaged in solution pillows. Requires CDTA Reagent, Cat. No. 2589636. 100 Tests.

Спецификации

Method Number 1 Chemistry:	8374
Бр. тестове:	100
Метод, име:	Chlorophosphonazo
Обем / опаковка, размер:	100 /pk
Параметър:	Calcium and Magnesium Hardness