



Be Right™



## EDTA Digital Titrator Cartridge, 0.1428 M

Продукт №:

1496001

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Ready-to-use titrant solution for the analysis of Total Hardness, 1-16 G.d.H. CaCO<sub>3</sub>.

For the determination of Calcium and Total Hardness. 0.1428 M. Digits x 0.1785 = mg/L CaCO<sub>3</sub> in 100 mL, Digits x 0.01 = Gdh.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Спецификации

Бр. тестове:	40 - 100
Измерване, диапазон:	1 - 16 G,d,h,
Метод, име:	EDTA (CalVer® 2)
Параметър:	Hardness, Calcium
Платформа:	Цифров титратор
Разтвор за титруване:	0.1428 M EDTA