



Tensette Pipette 1.0 - 10.0 mL

Продукт №:

1970010

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

Perform routine laboratory pipetting quickly, economically, and accurately with Hach's automatic Tensette Pipette. Two models dispense from 0.1 - 1 mL or from 1 - 10 mL in 0.1 or 1.0 mL increments, respectively. The pipette saves time, reduces glassware cleanup and assures accuracy. A Set includes one pipette, a manual and 50 disposable tips.

Спецификации

Обем / опаковка, размер:

1 each