



**Be Right™**



## **REF601 Reference electrode, Hg/Hg<sub>2</sub>SO<sub>4</sub>**

**Продукт №:**

**E21M012**

**ОСТАРЯЛ АРТИКУЛ**

Този артикул вече не е наличен.

Предложени артикули за замяна:

- E21M009
- B40A550

Reference electrode for chloride-free samples. With banana plug. Liquid junction: porous pin, salt bridge solution: sat. K<sub>2</sub>SO<sub>4</sub>. Temperature range 0-60°C.