



Hardness (Calcium) Reagent Set, LR, Digital Titrator

Продукт №:

2447200

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Reagent set for the determination of Calcium Hardness by EDTA Digital Titrator Titration Method (Hach method 8204). Range: 10-160 mg/L as CaCO₃. Set contains CalVer 2 Calcium Indicator Powder Pillows, Potassium Hydroxide Solution (8 N) and 0.0800 M (EDTA) Digital Titrator Cartridge for approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Спецификации

Бр. тестове:	40-100
Измерване, диапазон:	10 - 160 mg/L CaCO ₃
Метод, име:	EDTA Titration
Метод, номер:	8204
Параметър:	Hardness (Calcium)
Платформа:	Цифров титратор
Специална функция:	Reagent set, LR, Digital Titrator