



Натриев хлорид, 10 000 $\mu\text{S}/\text{cm}$, 500 mL

Продукт №: 2972249
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 10000 $\mu\text{S}/\text{cm}$ ± 100 $\mu\text{S}/\text{cm}$ at 25°C, 5.556 g/L ± 0.03 g/L as Sodium Chloride (NaCl), (1 $\mu\text{S}/\text{cm}$ \sim 0.5 mg/L TDS), 500 mL bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Спецификации

Package settings:	Bottle
Концентрация:	5,556 g/L $\pm 0,03$ g/L NaCl
Обем / опаковка, размер:	500 mL
Описание:	Натриев хлорид, 10 000 $\mu\text{S}/\text{cm}$, 500 mL
Параметър:	Conductivity
Платформа:	Разтвор
Условия на съхранение:	10 - 25 °C (protect from light)