



Radiometer Analytical XR440 Reference Electrode - New Hg free RedRod Reference with fiber Junction, Ground Joint, Screw Cap

Продукт №:

B20B440

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

The Radiometer XR440 Reference Electrode has a glass body, Red Rod reference system, refillable electrolyte (KNO_3 , 1 M), ground joint, and fiber junction. The XR440 reference probe is recommended for Voltammetry applications. The XR440 is a screw cap electrode and requires an external connection cable for connection to a meter.

The right probe for your application.

To meet your analysis needs, Radiometer Analytical reference electrodes are available with a variety of reference systems and liquid junctions. There is a Radiometer electrode for every application.

Спецификации

Гаранция:

12 месеца