



HQ14D Digital Conductivity meter kit, Cond. cell, Std., 1m

Продукт №:

HQ14D.99.201000

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Rugged portable meters for use in the field and plant.

Portable Conductivity meter. Basic instrument with Standard Conductivity cell, cable length 1m.

The digital meter/electrode system combines reliability, flexibility and ease of use. Interchangeable Intellical electrodes are recognised automatically and store all relevant data. Rugged, nearly indestructible outdoor versions with multiple cable lengths enable measurements even at previously inaccessible locations.

Replaced by NEW HQ Series.

Outstanding reliability and simple handling

Versatile electrodes for all applications, e.g. wastewater, drinking water, process water

Reliable results from inaccessible measurement locations and over long distances

Full GLP data management

Competent support through innovative technology and comprehensive applications know-how

Спецификации

AC, USB операции:	Optional
GLP функции:	Date; Time; Sample ID; Operator ID
IP клас на корпуса:	IP67
Parameters:	Conductivity
	NA
	NA
PC data transfer software :	HQD Series Meter Data Transfer Utility
Printer:	Yes, Optional
Salinity resolution:	0.01 (ppt) (‰)
Sensor A:	CDC40101 (1 m)
Sensors:	Sensor A: CDC40101
	Sensor B: NA

	Sensor C:
TDS измерване:	0 - 50,0 mg/L
TDS резолюция:	0.01 mg/L up - 0.1 g/L upon measuring range
Type of measurement:	Field Measurements
Вклчени сонди:	CDC40101
Входове:	1
Гаранция:	36 месеца
Данни, експортиране:	Download via USB connection to PC or flash memory device. Automatically transfer entire data log or as readings are taken. Power supply adapter required.
Данни, съхранение:	Automatic, GLP ISO compliant reading data stored with batteriescalibration details.
Данни, съхранение / памет:	500 records/FIFO
Дисплей, тип:	240 x 160 pixel LCD with batteriesbacklight illumination
Дисплей, функция за заключване:	Continuous / Auto-stabilization ("press to read") / At Interval
Дисплей: течнокристален (фоново осветление):	Detailed mode/Large mode
Езици:	English, German, French, Italian, Spanish, Danish, Dutch, Polish, Portuguese, Turkish, Finnish, Czech, Russian.
Измерване метод:	Probe specific programmed method settings
Измерване на соленост:	0 - 42 (ppt) (‰)
Инструмент:	Portable
Интервали / предупреждения / напомняния при калибриране:	Изключено, избираем от 2 часа до 7 дни
Интерфейс:	Soft Touch Keypad
Какво има в кутията:	Meter package includes HQ14D Portable Meter; 4 AA batteries; CDC401 Laboratory Conductivity Cell, 1 m Cable;
	Sodium Chloride Standard Solution, 491 mg/L NaCl (1000 µS/cm), 100 mL; quick-start guide and user manual.
Калибриране на електрод CDC:	Demal (1D/ 0.1D/ 0.01D);
	Molar (0.1M/ 0.01M/0.001M);
	NaCl (0.05%; 25µS/cm; 1000µS/cm; 18mS/cm);
	Standard sea water;
	User defined
Комплект?:	Да
Модел:	HQ14D – EC/1 канал
Обем на доставката:	Meter package includes HQ14D Portable Meter; 4 AA batteries; CDC401 Laboratory Conductivity Cell, 1 m Cable;
	Sodium Chloride Standard Solution, 491 mg/L NaCl (1000 µS/cm), 100 mL; quick-start guide and user manual.
Околна среда, условия: относителна влажност:	90 % relative humidity (non-condensing)
Околна среда, условия: температура:	0 - 60 °C (32 - 140 °F)
Параметър:	Conductivity/Total Dissolved Solid (TDS)/Salinity/Resistivity
Проводимост, измерване: корекция на температурата:	None; Linear; NaCl Non-Linear Natural Water.
Проводимост: измерване:	0.01 µS/cm - 200.0 mS/cm

Проводимост: измерване при стабилни показания:	Yes
Проводимост: резолюция:	0.01 µS/cm - 0.1 mS/cm upon selected measuring range
Проводимост: точност:	± 0.5 % in (1µS/cm -200mS/cm) range
Сертификати за съответствие:	CE.WEEE
Стандартни разтвори за специално калибриране:	Yes
Съдържание:	No
съобщения за грешки при работа:	Clear text error messages displayed
Тегло:	342 g (0.75 lb) without batteries; 430 g (0.95 lb) with batteries
Температура, измерване:	°C or °F
Температура, резолюция:	0,1 °C
Температура, точност:	± 0.3 °C
Тип на сондата:	Intellical Digital SMART Standard Laboratory or Rugged Field

Какво има в кутията

Meter package includes HQ14D Portable Meter; 4 AA batteries; CDC401 Laboratory Conductivity Cell, 1 m Cable; Sodium Chloride Standard Solution, 491 mg/L NaCl (1000 µS/cm), 100 mL; quick-start guide and user manual.