



## Beaker, polypropylene, low form, 50 mL

Продукт №:

108041

BGN Цена:

Свържете се с нас

Налично

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Спецификации

Диаметър, външен:	51.1 mm with spout
Единица:	1 each
Капацитет:	50 mL
Материал:	Polypropylene, Low Form
Материал на тялото (измервателния уред):	Polypropylene
Обем / опаковка, размер:	1 each