



Beaker, polypropylene, low form, 1000 mL

Продукт №:	108053
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Спецификации

Диаметър, външен:	129.8 mm with spout
Единица:	1 each
Капацитет:	1000 mL
Материал:	Polypropylene, Low Form
Материал на тялото (измервателния уред):	Polypropylene