



## Beaker, polypropylene, low form, 100 mL, 12/pk

Продукт №:

108072

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Спецификации

Диаметър, външен:	63.0 mm with spout
Единица:	12 /pk
Капацитет:	100 mL
Материал:	Polypropylene, Low Form
Материал на тялото (измервателния уред):	Polypropylene