



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Продукт №:

108076

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Спецификации

Диаметър, външен:	83.6 mm with spout
Единица:	6 /pk
Капацитет:	250 mL
Материал:	Polypropylene, Low Form
Материал на тялото (измервателния уред):	Polypropylene