



## Flask, volumetric, polypropylene, 100 mL

|                                |                   |
|--------------------------------|-------------------|
| Продукт №:                     | 1406042           |
| BGN Цена:                      | Свържете се с нас |
| Не е показана дата на доставка |                   |

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

---

### Спецификации

|                                          |                |
|------------------------------------------|----------------|
| Единица:                                 | 1 each         |
| Капацитет:                               | 100 mL         |
| Материал на тялото (измервателния уред): | Polypropylene  |
| Сертификати:                             | Not Applicable |