



Flask, volumetric, polypropylene, 250 mL

Продукт №:

1406046

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Спецификации

Единица:	1 each
Капацитет:	250 mL
Материал на тялото (измервателния уред):	Polypropylene
Сертификати:	Not Applicable