



Hardness, iron, and pH test kit, model HA-62B

Продукт №:

183702

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Model HA-62B

Includes color disc and color comparator for colorimetric comparison of iron and pH, as well as a drop-count titrator for hardness determination. Test kit produces a wide range of pH tests.

Combines frequently measured parameters in one rugged kit

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Бр. тестове:	100
Модел:	HA-62B
Параметър:	Hardness, Iron, and pH (4-10)
Платформа:	Цветен диск
Приложение:	Water Quality