



Cyanuric acid test kit, model CY-3

Продукт №: 185102
BGN Цена: Свържете се с нас
Не е показана дата на доставка

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Smallest increment:	20
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	20 - 100 mg/L
Метод, име:	Turbidimetry
Модел:	CY-3
Модел на тестов комплект:	CY-3
Параметър:	Cyanuric Acid
Приложение:	Pool and Spa