



Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Продукт №:	190902
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

Спецификации

Body Material:	Borosilicate Glass
Package size UOM:	1pk
Капацитет:	60 mL
Материал:	Borosilicate Glass
Обем / опаковка, размер:	1 each