



Test kit tannin/lignin

Продукт №:

193701

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Опасен



A product of decomposed plant material, tannin occurs in natural waters and is used for boiler and cooling tower water treatment, dyeing, ink manufacture, tanning, and paper sizing. Lignin, also a plant product, is found in paper pulp manufacturing waste. Tannin and lignin are analyzed together and reported as mg/L tannic acid. This kit utilizes the simple, low-cost color disc method - the most accurate visual test kit method.

Tannin and lignin are analyzed together and reported as mg/L tannic acid.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Parameter:	Tannin/Lignin
Smallest increment:	0.5
	5
Бр. тестове:	100
Вид на корпуса:	R
Измерване, диапазон:	0 - 15 mg/L as Tannic Acid
Метод, име:	Color disc/Tyrosine
Модел:	Tannin/Lignin
Модел на тестов комплект:	TA-3
Параметър:	Tannin/Lignin - As tannic acid
Платформа:	Цветен диск
Стандартен обхват на измерване:	0 - 15 mg/L as Tannic Acid; 0 - 150 mg/L