



## Bag, sterile, WhirlPak w/ dechlorinating agent, 100 mL, 25/pk

Продукт №: 2075325  
BGN Цена: Свържете се с нас  
Не е показана дата на доставка

Specially designed for chlorinated water samples (Complies with USEPA and Standard Methods criteria). Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Contains non-nutritive pill with 10 mg sodium thiosulfate (will dechlorinate up to 10 mg/ $\text{Cl}_2$  in 100 mL sample) Four ounce and 100 mL fill-lines. "Write-on" strip accepts waterproof ink.

---

### Спецификации

Капацитет:	177 mL
Обем / опаковка, размер:	25