



Test kit, chloride, model CD-51

Продукт №:

208601

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Опасен



Found in nearly all natural waters, chloride affects human taste above 250 mg/L. High levels inhibit plant growth, and many industrial processes need to limit chloride concentration. The model CD-DT, multi-range kit, is supplied with the demineralizer bottle for high-range sample dilution. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Smallest increment:	500
	5000
Бележка под линия:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Бр. тестове:	100
Вид на корпуса:	B
Измерване, диапазон:	500 - 10000 mg/L & 5000 - 100000 mg/L Cl ⁻
Метод, име:	Drop count titration/
	Silver nitrate
Модел на тестов комплект:	CD-51
Параметър:	Chloride - As Cl ⁻
Платформа:	Drop Count Titration