



Bottle, Sampling, Low Density, Polyethylene, 250 mL, 12/pk

Продукт №:	2087076
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Спецификации

Body Material:	Low Density Plastic
Package size UOM:	12pk
Капацитет:	250 mL
Материал:	LDPE
Обем / опаковка, размер:	12 /pk