



Bottle, Sampling, Low Density, Polyethylene, 1 L, 6/pk

Продукт №:	2087083
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Спецификации

Body Material:	Low Density Plastic
Package size UOM:	6pk
Капацитет:	1000 mL
Материал:	LDPE
Обем / опаковка, размер:	6 /pk