



**Be Right™**



## Rack, drying

Продукт №:

**2089000**

EUR Цена:

Свържете се с нас

Не е показана дата на доставка

Epoxy-coated, steel-wire drying rack for test tubes, vials, flasks, beakers, and other glassware. Three rows of hair-pin looped supports to hold large-diameter necks. Fourteen loops total. Top row has 12 straight-pin for smaller tubes and glassware. Polystyrene tray fits into bottom rails to catch draining water. Dimensions: 43 x 17 x 48 cm (17 x 6.7 x 19"). Free-standing or wall mount option.

---

### Спецификации

Дълбочина:	48 cm
Дължина:	430 mm
Обем / опаковка, размер:	1 each