



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 250 mL

Продукт №:

2089676

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Спецификации

Единица:	6/pk
Капацитет:	250 mL
Материал:	Glass Screw Cap