



Flask, Erlenmeyer, polymethylpentene capacity 250 mL, 4/pk

Продукт №: 2089876
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Спецификации

Вид:	Erlenmeyer
Единица:	4/pk
Капацитет:	250 mL
Материал:	Polymethylpentene