



Weight Set, Class M-2

Продукт №: 2101100
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Standard weight set. 16-piece, with stainless steel forceps. Eight fractional weights: 10, 20, 50, 100, 200, 200, 200, and 500 mg. (Two fractional weights, 200 mg and 500

mg, are stainless steel, the rest are aluminum) Eight brass gram weights: 1, 2, 2, 5, 10, 10, 20, and 50 grams. Polypropylene case.

Спецификации

Какво има в кутията:	Includes: Brass weights and polypropylene case.
Обем / опаковка, размер:	1 each

Какво има в кутията

Includes: Brass weights and polypropylene case.