



Be Right™



Centrifuge Tube, Conical Bottom, 10/pk

Продукт №:

2115839

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Centrifuge tube with conical bottom. Polycarbonate with molded-in 0.1mL graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is also used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Autoclavable. Nalge.

Спецификации

Дължина:	120 mm
Капацитет:	15 mL
Обем / опаковка, размер:	10 /pk