



Cylinder, graduated, polymethylpentene, 250 mL $\pm$ 1.4 mL, 2.0 mL divisions

Продукт №: 217246
BGN Цена: Свържете се с нас
Не е показана дата на доставка

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated "To Contain/To Deliver" at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Спецификации

Вид:	Две метрични стойности
Измерване метод:	To Deliver/To Contain
Капацитет:	250 mL
Материал на тялото (измервателния уред):	Polymethylpentene
Обем / опаковка, размер:	1 each
Сертификати:	Class B