



Dish, Petri, 15 x 100 mm

Продукт №:	2178900
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Widely used in microbiology, these standard-size polystyrene petri dishes are presterilized, stackable, and disposable. Dishes are packaged 24 per sleeve, with a case containing 25 sleeves (600 dishes total).

Standard-size polystyrene petri dishes

Presterilized, stackable, and disposable

600 dishes total

Спецификации

Обем / опаковка, размер:	600
--------------------------	-----