



## Test kit, formaldehyde, model FM-1

Продукт №:

2183100

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

**Utilises drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Спецификации

|                           |                                          |
|---------------------------|------------------------------------------|
| Smallest increment:       | 0.05%                                    |
|                           | 0.5%                                     |
| Бр. тестове:              | 100                                      |
| Вид на корпуса:           | D                                        |
| Измерване, диапазон:      | 0,05 - 1% & 0,5 - 10%                    |
| Метод, име:               | Drop count titration/<br>Thymolphthalein |
| Модел на тестов комплект: | FM-1                                     |
| Параметър:                | Formaldehyde - As $\text{CH}_2\text{O}$  |
| Платформа:                | Drop Count Titration                     |