



Test kit, general boiler water, drop count titration, model AL-94 MG-L

Продукт №:

222403

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Опасен



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus to perform analysis in the field

Uses drop count titration and color disc methods

Спецификации

Модел:	AL-94 MG-L
Параметър:	Drop Count Titration
Платформа:	Цветен диск
Приложение:	Boiler Water