



## Test kit, nitrogen , model NI-8

Продукт №:

224100

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Опасен



Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

**One of Hach's most popular test kits**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Спецификации

|                           |                                               |
|---------------------------|-----------------------------------------------|
| Smallest increment:       | 0.1                                           |
| Бр. тестове:              | 100                                           |
| Вид на корпуса:           | B                                             |
| Измерване, диапазон:      | 0 - 3,0 mg/L NH <sub>3</sub> -N               |
| Метод, име:               | Color disc/Nessler reagent                    |
| Модел на тестов комплект: | NI-8                                          |
| Параметър:                | Nitrogen, Ammonia, mid range - For freshwater |
| Платформа:                | Цветен диск                                   |