



Test kit, chlorine (total), model CN-65

Продукт №:

225401

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Smallest increment:	0.2
	1
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	0,2 - 4,0 mg/L & 1 - 20 mg/L
Метод, име:	Drop count titration/
	Thiosulfate
Модел на тестов комплект:	CN-65
Най-малки инкременти, стъпки:	1
Параметър:	Chlorine, mid and high range - As Total Cl ₂
Платформа:	Drop Count Titration