



SA filter, cone, diameter 110 mm, soil analysis grade, 500/pk

Продукт №: 2263356
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Designed for sandy soils, SA filters are conically folded and have been optimized for use with the common extraction solutions used in soil nutrient analysis (Mehlich 1, Mehlich 3, Morgan, Modified Morgan, Bray, and Olsen). SA Filters provide the best combination of retention rating, filtration speed and purity requirements for soil nutrient testing (P, K, Ca, Mg). Medium filter retention.

Спецификации

Диаметър:	110 mm
Материал:	Cellulose
Обем / опаковка, размер:	50 /pk