



Filter, acrocap, pore size 0.2 µm

Продукт №: 2323810
BGN Цена: Свържете се с нас
Не е показана дата на доставка

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

Спецификации

Filter Type:	Filter, acrocap
Packing:	NA
Диаметър:	NA
Пори, размер:	0.2 µm