



Be Right™



Bottle, sampling, polycarbonate, 500 mL, 4/pk

Продукт №:	2324539
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Transparent polycarbonate bottle with cap features wide mouth (3.8cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Спецификации

Body Material:	Polycarbonate
Капацитет:	500 mL
Обем / опаковка, размер:	4 /pk