



Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

Продукт №:

2349500

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Бр. тестове:	50
Вид на корпуса:	R
Измерване, диапазон:	0 - 300 µg/L DEHA (and other Oxygen Scavengers)
Метод, име:	Color disc/Iron reduction
Модел:	OS-1
Модел на тестов комплект:	OS-1
Параметър:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Платформа:	Цветен диск